



## Kit Quick-TB

**Code product:** 361560-0000

### Composition:

- Carbolic Fuchsin RAL
- Armand Solution

Code: 365240-0125

1 x 125 ml

Code: 3601002A125

1 x 125 ml

### Applications:

#### Bacteriology:

Staining with Kit Quick-TB allows to detect mycobacteria, also called A.F.B. (Acid-Alcohol-Fast Bacilli) through a cold and fast staining technique.

Carbolic Fuchsin RAL stains mycobacteria pink.

Armand Solution stains the background of the preparation blue.

Updated: 25/04/2013