



Kit Fluo-RAL

Code product : 359000-0000

Composition:

- | | |
|---|------------------------------|
| • Fixative Fluo-RAL | code: 362115-0200 1 x 200 mL |
| • Auramine Fluo-RAL | code: 362325-0100 1 x 100 mL |
| • Carbolic Solution for Auramine Fluo-RAL | code: 362530-0100 1 x 100 mL |
| • Degommier Discolouring Fluo-RAL | code: 362210-0200 1 x 200 mL |
| • Thiazine Red Fluo-RAL | code: 362780-0100 1 x 100 mL |
| • Carbolic Solution for Thiazine Red Fluo-RAL | code: 362610-0100 1 x 100 mL |

Applications:

Bacteriology:

Staining with Kit Fluo-RAL allows to detect mycobacteria, also called A.F.B. (Acid-Alcohol-Fast Bacilli) through a technique by fluorescence according to Degommier. Auramine Fluo-RAL, extemporaneously mixed with Carbolic Solution for Auramine Fluo-RAL, stains mycobacteria fluorescent greenish yellow. Degommier Discolouring Fluo-RAL discolour all bacteria except mycobacteria. Thiazine Red Fluo-RAL, extemporaneously mixed with Carbolic Solution for Thiazine Red Fluo-RAL, stains the background of the preparation red.

Updated: 25/04/2013